

Product datasheet for **MR208523**

Btbd3 (BC018174) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Btbd3 (BC018174) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Btbd3
Synonyms:	MGC25591, mKIAA0952
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR208523 representing BC018174
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGTAGATGATAAAGAAAAGAACATGAAATGTCTCACCTTCTTCTTGATGCTTCCAGAGACAGTAAAGA
 ACAGGTCCAAAAAGGTTCAAAGAAAGCCAACAGCAGCGGTGGCGCGCGCGCGGCAGCGTTGGCAG
 TGGGAGCAGCAAGTTGCCCCAGTTTGTATGAAATAATTACCTTGAAGACTAAAAAGAAGAAGATG
 GCTGCTGATATATTCCCCGTAAGAAGCCAGCCAACTCCAGCAGCACCCTGTCCAGCAGCAGCACCAGC
 ACCAGCACAACCTCTGTAACAACAACCTCATCCCGCCCCCAACTGGCAGGGTCTTTATCCACCATCCG
 AGAAAGAAATGCGGTGATGTTCAATAATGACTTGATGGCAGACGTACATTTTGTGGTTGGGCCACCAGGT
 GGGACTCAGCGGTTGCCAGGACACAAATACGTCTCGTGTGGGAGCTCAGTGTTCATGCAATGTTTT
 ACGGAGAACTCGCTGAGGACAAAGATGAAATCCGAATCCAGACGTGAACTGTGCCTTCTCGCCAT
 GCTGAAATACATCTATTGTGATGAAATTGACTTGGCTGCAGACACAGTGTGGCCACACTGTATGCTGCC
 AAAAGTACATAGTTCTCACCTCGCCAGAGCCTGTGTTAATTTCTGGAGACCAGCCTAAGCGCCAAGA
 ACGCCTGCGTACTCTCTCCAGAGCTGCCTGTTTGAAGAGCCGACCTGACACAGCGTTGCTGGGAGGT
 GATTGATGCCCAGGCCGAGTTAGCGCTCAAGTCTGAGGGATTCTGTGACATCGACTTCCAGACACTCGAA
 AGTATCTCCGCAAGGAAACGCTGAATGCCAAAGAAATAGTGGTTTTGAGGCAGCCCTGAAGTGGGCTG
 AGGTAGAATGCCAGCGCAAGATCTGGCCCTGAGCATTGAAAAATAACGCAAGGTCTAGGAAAGGCCCT
 GTACTTGATCCGCATCCCCACCATGGCCCTGGATGACTTGGCAATGGTGTGCACAGTCTGGGGTATTG
 ACTCTCAATGAGACCAATGACATCTTCTCTGGTACACTGCATCCAAAAACCAGAGCTGCAGTTTGTGA
 GTAAAGCCCGCAAGGGCTGGTCCCCAGCGTTGCCACCGATTTCAGTCATGTGCCTATCGGAGCAACCA
 GTGGCGCTATCGGGTTCGCTGTGACAGCATCCAGTTTGTGTTGATAAGAGGGTGTTCATCGCAGGCTTT
 GGGCTGTACGGCTCCAGCTGTGGCTCTGCCGAGTACAGCGCAAAGATCGAACTCAAACGGCAGGGCGTTG
 TCCTGGGCCAGAACTTGAGCAAGTACTTCTCAGACGGATCCAGCAATACCTTCCCTGTGTGGTTTGA
 TCCAGTGCAGATTGAGCCGGACACCTTCTACACAGCCAGTGTAGTGTGGACGGCAATGAAGTACAGTAC
 TTTGGGCAGGAAGGCATGACGGAAGTTCAGTGTGGCAAGGTGACTGTCCAGTTCAGTGTCTCTCAGACA
 GCACCAATGGCACTGGTGTACAGGGAGGCAGATCCCCGAATCATATTCTACGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR208523 representing BC018174
 Red=Cloning site Green=Tags(s)

MVDDKEKNMKCLTFFLMLPETVKNRSKKGSKKANSSGGGGGGSVSGSSKLPPVCYEIITLKTKKKKKM
 AADIFPRKKPANSSSTTVQQHQHQHNLNNLIPAPNWQGLYPTIRERNAVNFNDLMADVHFVVGPPG
 GTQRLPGHKYVLAVGSSVFHAMFYGELAEDKDEIRIPDVEPAFLAMLKYIYCDEIDLAADTVLATLYAA
 KKYIIVPHLARACVNFLETSLSAKNACVLLSQSCLFEEDLTQRCWEVIDAQELALKSEGFCIDFQTL
 SILRRETLNAKEIVVFEALNWAEEVCQRQDLALSIENKRKVLGKALYLIRIPTMALDDFANGAAQSGVL
 TLNETNDIFLWYASKKPELQFVSKARKGLVPQRCHRFQSCAYRSNQWRYRGRCDSIQFAVDKRVFIAGF
 GLYSSCGSAEYSAKIELKRQGVVLGQNLISKYSDGSSNTFPVWFEYPVQIEPDTFYTASVVLGDNELSY
 FGQEGMTEVQCGKVTQVQCSSDSTNGTGVQGGQIPELIFYA

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: BC018174

ORF Size: 1596 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC018174.1](#)


RefSeq Size: 2436 bp

RefSeq ORF: 1598 bp

Locus ID: 228662

Cytogenetics:	2 F3
MW:	89.3 kDa
Gene Summary:	Acts as a key regulator of dendritic field orientation during development of sensory cortex. Also directs dendrites toward active axon terminals when ectopically expressed. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR208523