

Product datasheet for **MR208530L4V**

Eya2 (NM_010165) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Eya2 (NM_010165) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Eya2
Synonyms:	Eab1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010165
ORF Size:	1599 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208530).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010165.2
RefSeq Size:	2460 bp
RefSeq ORF:	1599 bp
Locus ID:	14049
UniProt ID:	O08575
Cytogenetics:	2 85.83 cM
Gene Summary:	This gene encodes a member of the eyes absent protein family. The encoded protein is a tyrosine phosphatase which acts as a transcriptional activator during development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2012]



[View online »](#)