

## Product datasheet for **MR208549L3V**

### Neto1 (NM\_144946) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Neto1 (NM_144946) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Neto1
Synonyms:	AI851453; Btcl1; C130005O10Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144946
ORF Size:	1599 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208549).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_144946.4</a>
RefSeq Size:	3531 bp
RefSeq ORF:	1602 bp
Locus ID:	246317
UniProt ID:	<a href="#">Q8R4I7</a>
Cytogenetics:	18 E4



[View online »](#)

**Gene Summary:**

Involved in the development and/or maintenance of neuronal circuitry. Accessory subunit of the neuronal N-methyl-D-aspartate receptor (NMDAR) critical for maintaining the abundance of GRIN2A-containing NMDARs in the postsynaptic density. Regulates long-term NMDA receptor-dependent synaptic plasticity and cognition, at least in the context of spatial learning and memory.[UniProtKB/Swiss-Prot Function]