

Product datasheet for **MR208594L3V**

Pcdhb22 (BC094239) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Symbol: Pcdhb22

Synonyms: PcdhbV, Pcdhb15

Mammalian Cell Selection: Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK

ACCN: BC094239

ORF Size: 1611 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR208594).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [BC094239](#), [AAH94239](#)

RefSeq Size: 2287 bp

RefSeq ORF: 1613 bp

Locus ID: 93893

Cytogenetics: 18 B3

