

Product datasheet for **MR208611L4V**

Cspg5 (NM_013884) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cspg5 (NM_013884) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cspg5
Synonyms:	CAL; Caleb; NG; Ngc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013884
ORF Size:	1620 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208611).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013884.3
RefSeq Size:	3800 bp
RefSeq ORF:	1620 bp
Locus ID:	29873
UniProt ID:	Q71M36
Cytogenetics:	9 F2



[View online »](#)

Gene Summary:

This gene encodes a chondroitin sulfate proteoglycan. The encoded protein has been termed a 'part-time' proteoglycan, as chondroitin sulfate chains appear to be attached to the protein in the developing but not the adult cerebellum and retina. It is thought that this protein plays roles in dendrite branching and synapse formation. [provided by RefSeq, Oct 2009]