

Product datasheet for **MR208638L3V**

Net1 (NM_001047159) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Net1 (NM_001047159) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Net1
Synonyms:	mNET1; Net1a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001047159
ORF Size:	1626 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208638).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001047159.1 , NP_001040624.1
RefSeq Size:	3874 bp
RefSeq ORF:	1626 bp
Locus ID:	56349
UniProt ID:	Q9Z206
Cytogenetics:	13 A1
Gene Summary:	Acts as guanine nucleotide exchange factor (GEF) for RhoA GTPase. May be involved in activation of the SAPK/JNK pathway. Stimulates genotoxic stress-induced RHOB activity in breast cancer cells leading to their cell death.[UniProtKB/Swiss-Prot Function]



[View online »](#)