

Product datasheet for **MR208653L3V**

Nsrp1 (NM_001012309) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nsrp1 (NM_001012309) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nsrp1
Synonyms:	8030451K07; AI851076; AL022898; Ccdc55; NSpr70; NSrp70
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001012309
ORF Size:	1626 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208653).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001012309.2 , NP_001012309.1
RefSeq Size:	3132 bp
RefSeq ORF:	1629 bp
Locus ID:	237859
UniProt ID:	Q5NCR9
Cytogenetics:	11 B5
Gene Summary:	RNA-binding protein that mediates pre-mRNA alternative splicing regulation. [UniProtKB/Swiss-Prot Function]



[View online »](#)