

Product datasheet for **MR208685L4V**

Leprel2 (BC003726) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Leprel2 (BC003726) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Leprel2
Synonyms:	BC016431; Grcb; Leprel2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC003726
ORF Size:	1635 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208685).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC003726 , AAH03726
RefSeq Size:	2289 bp
RefSeq ORF:	1637 bp
Locus ID:	14789
Cytogenetics:	6 59.17 cM



[View online »](#)

Gene Summary:

Part of a complex composed of PLOD1, P3H3 and P3H4 that catalyzes hydroxylation of lysine residues in collagen alpha chains and is required for normal assembly and cross-linking of collagen fibrils (PubMed:27119146). Required for normal hydroxylation of lysine residues in type I collagen chains in skin, bone, tendon, aorta and cornea (PubMed:28115524). Required for normal skin stability via its role in hydroxylation of lysine residues in collagen alpha chains and in collagen fibril assembly (PubMed:27119146, PubMed:28115524). Apparently not required for normal prolyl 3-hydroxylation on collagen chains, possibly because it functions redundantly with other prolyl 3-hydroxylases (PubMed:28115524).[UniProtKB/Swiss-Prot Function]