

Product datasheet for **MR208695L4V**

Ddx56 (NM_026538) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ddx56 (NM_026538) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ddx56
Synonyms:	2600001H07Rik; D11Ert619e; Noh61
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026538
ORF Size:	1641 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208695).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026538.3 , NP_080814.1
RefSeq Size:	3287 bp
RefSeq ORF:	1641 bp
Locus ID:	52513
UniProt ID:	Q9D0R4
Cytogenetics:	11 3.94 cM
Gene Summary:	May play a role in later stages of the processing of the pre-ribosomal particles leading to mature 60S ribosomal subunits. Has intrinsic ATPase activity (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)