

Product datasheet for **MR208698L3V**

Odf2l (NM_001162539) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Odf2l (NM_001162539) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Odf2l
Synonyms:	4733401D09Rik; 9630045K08Rik; D3Ert250e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001162539
ORF Size:	1641 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208698).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001162539.1 , NP_001156011.1
RefSeq Size:	2296 bp
RefSeq ORF:	1641 bp
Locus ID:	52184
UniProt ID:	Q9D478
Cytogenetics:	3 69.36 cM
Gene Summary:	Acts as a suppressor of ciliogenesis, specifically, the initiation of ciliogenesis. [UniProtKB/Swiss-Prot Function]


[View online »](#)