

## Product datasheet for **MR208739L3V**

### Aspscr1 (NM\_026877) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Aspscr1 (NM_026877) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Aspscr1
Synonyms:	1190006K01Rik; ASPC; ASPCR1; ASPL; ASPS; RCC17; TUG
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026877
ORF Size:	1653 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208739).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026877.1</a>
RefSeq Size:	1872 bp
RefSeq ORF:	1653 bp
Locus ID:	68938
UniProt ID:	<a href="#">Q8VBT9</a>
Cytogenetics:	11 E2
Gene Summary:	Enhances VCP methylation catalyzed by VCPKMT (By similarity). Tethering protein that sequesters GLUT4-containing vesicles in the cytoplasm in the absence of insulin. Modulates the amount of GLUT4 that is available at the cell surface.[UniProtKB/Swiss-Prot Function]



[View online »](#)