

Product datasheet for **MR208777L4V**

Mettl16 (NM_026197) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mettl16 (NM_026197) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mettl16
Synonyms:	2610100D03Rik; 2810013M15Rik; A830095F14Rik; AI846133; Mett10d
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026197
ORF Size:	1662 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208777).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026197.1
RefSeq Size:	2718 bp
RefSeq ORF:	1662 bp
Locus ID:	67493
UniProt ID:	Q9CQG2
Cytogenetics:	11 B5



[View online »](#)

Gene Summary:

RNA N6-methyltransferase that methylates adenosine residues at the N(6) position of a subset of RNAs and is involved in S-adenosyl-L-methionine homeostasis by regulating expression of MAT2A transcripts (PubMed:29262316, PubMed:30197299). Able to N6-methylate a subset of mRNAs and U6 small nuclear RNAs (U6 snRNAs) (By similarity). In contrast to the METTL3-METTL14 heterodimer, only able to methylate a limited number of RNAs: requires both a 5'UACAGAGAA-3' nonamer sequence and a specific RNA structure (By similarity). Plays a key role in S-adenosyl-L-methionine homeostasis by mediating N6-methylation of MAT2A mRNAs, altering splicing and/or stability of MAT2A transcripts: in presence of S-adenosyl-L-methionine, binds the 3' UTR region of MAT2A mRNA and specifically N6-methylates the first hairpin of MAT2A mRNA, impairing MAT2A expression (PubMed:29262316, PubMed:30197299). In S-adenosyl-L-methionine-limiting conditions, binds the 3' UTR region of MAT2A mRNA but stalls due to the lack of a methyl donor, preventing N6-methylation and promoting expression of MAT2A (PubMed:29262316). In addition to mRNAs, also able to mediate N6-methylation of U6 small nuclear RNA (U6 snRNA): specifically N6-methylates adenine in position 43 of U6 snRNAs (By similarity). Also able to bind various lncRNAs (By similarity). Specifically binds the 3'-end of the MALAT1 long non-coding RNA (By similarity) (PubMed:29262316, PubMed:30197299).[UniProtKB/Swiss-Prot Function]