

Product datasheet for **MR208781L3V**

Zcchc5 (Rtl3) (NM_199468) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Zcchc5 (Rtl3) (NM_199468) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rtl3
Synonyms:	D430021I08Rik; Gm375; Mar3; Mart3; ZHC5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_199468
ORF Size:	1659 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208781).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_199468.1
RefSeq Size:	3019 bp
RefSeq ORF:	1662 bp
Locus ID:	213436
UniProt ID:	Q6P1Y1
Cytogenetics:	X D
Gene Summary:	May function as a transcriptional regulator (Probable). Plays a role in postnatal myogenesis, may be involved in the regulation of satellite cells self-renewal (PubMed:27446912). [UniProtKB/Swiss-Prot Function]



[View online »](#)