

Product datasheet for **MR208796L4V**

Hnrnpl (NM_177301) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hnrnpl (NM_177301) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hnrnpl
Synonyms:	C79783; D830027H13Rik; Hnrpl
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_177301
ORF Size:	1758 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208796).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177301.5 , NP_796275.3
RefSeq Size:	2142 bp
RefSeq ORF:	1761 bp
Locus ID:	15388
UniProt ID:	Q8R081
Cytogenetics:	7 B1



[View online »](#)

Gene Summary:

Splicing factor binding to exonic or intronic sites and acting as either an activator or repressor of exon inclusion (PubMed:22523384). Exhibits a binding preference for CA-rich elements. Component of the heterogeneous nuclear ribonucleoprotein (hnRNP) complexes and associated with most nascent transcripts. Associates, together with APEX1, to the negative calcium responsive element (nCaRE) B2 of the APEX2 promoter.[UniProtKB/Swiss-Prot Function]