

Product datasheet for **MR208812L3V**

Prosapip1 (BC058280) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Prosapip1 (BC058280) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Prosapip1
Synonyms:	Prosapip1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC058280
ORF Size:	1668 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208812).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC058280 , AAH58280
RefSeq Size:	3832 bp
RefSeq ORF:	1670 bp
Locus ID:	241638
Cytogenetics:	2 F1
Gene Summary:	May be involved in promoting the maturation of dendritic spines, probably via regulating SIPA1L1 levels at the postsynaptic density of synapses.[UniProtKB/Swiss-Prot Function]


[View online »](#)