

Product datasheet for **MR208965**

Eepd1 (NM_026189) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Eepd1 (NM_026189) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Eepd1
Synonyms:	2310005P05Rik; AI481005
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR208965 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGGAGCACACTGGGCTGCCATCGCTCTATCCCCAGAGATCCCTCGGACCTGTACATAACCGCAAGT
 TCAGCGCTGCCTGCAACTTCAGCAACATCCTGGTGAATCAGGAGCGGCTAAATATCAACACTGCCACTGA
 GGAGGAAGTATGACTCTGCCTGGGGTGACACGTGCAGTGGCCCCGAAGCATCGTGAATACCGAGAATAT
 ATCGGGCGGCTTCAAAAAGGTGGAAGATCTGGCACTTGTCAAGTGGTGTAGGGGCCACCAAAGTGAACAGG
 TCAAGTTTGAGATCTGTGTGAGCAGCAAGGTAAGTCTGCACAGCACTCTCCAGTTCCCTGAGGCGGGA
 TCTTCTGGCTGAGCAGCAGCCTCATCACCTGACCACCACTGTGCCCTCACTCCGCGAGTCAATATCAAC
 ACAGCCACCCCTTGCTCAGCTCATGAGTGTCCGGGGCTCAGCGAGAAGATGGCCCTCAGCATCGTGGACT
 ACCGCAGGGAACACGGTCTTTTCGAAGTGTGAGGACTTGGTGAAGTGGATGGAATCAATGCAGCCTT
 CCTGGACAGAATAAGACACCAGGTGTTTCGAGAGAGGTCCAGGCCCATCAACCCACACCAATGGGGT
 CTTACCTTCACTGCCAAGCCTCATCCCAGTCCCACATCTCTGAGCCTGCAGAGTGAAGACTTGGACTCTGC
 CACCAGGGGGGCCACGAGATTATCTCCATGAGGCCCTCGGTGGAGGCCTTTGGAGGTATGAGGGATGG
 GCGACCTGTGTTAGGTTGGCCACGTGGAACCTACAGGGCTGCTCTGTGAGAAGGCCAATAACCGGGG
 GTTCGAGAGGTGGTGTGCATGACGCTTCTGAAAAACAGCATCAAAGTCTAGCTGTGCAAGAGCTGCTTG
 ACAAAGAGGCCCTGAAAAAGTCTGTACAGAGCTGAACAGCCGATCCTGCCCAACATCCGCAAGTGAA
 GGGTCCCAGGATGCTGGAGGTCCATTGTTGCTGAGAAGCCTTCCAATCAGTCCAGAAGGGCCCTGC
 TACTCTGGATTCTATGGGACACAGCAGCAACGTGGAGCTGAGAGACATCCCCGGGCGGAGAGCTCAC
 CCAGAACCGTCACTAAGGCTGTGGGCCCCAGCCCTTCTCACCAGATTCAAGGTGGGAAGTAAATGA
 CTTAACCTGGTTAACCTCCAGCTGACAGCCCTGGCCCTCCAGGAGCTGAGAATTCAAGCAAGAAATCAC
 AGTGATGGCCACCGATTATTGAAGTTTGCATAACCTCCAGGAAACCTGAAAGGAGAAAAAGGATGTTG
 TTATTTTAGGGGATTTGGCCAAGGGCCAGATAGCAATGACTACGATATCCTAAGGAGAGAAAAATTTCA
 CCACCTAGTCCCGGCTCACACTTTCACCAACATCAGCACCAGGAATCCTCAAGGCTCAAGTCTGTGGAC
 AACATCTGGATCAGCAAAAGCTTAAAGAAAGTGTTCACAGGCTACTGGGCTGTTGTGAGAGAAGGGCTCA
 CAAACCTTGGATTCCAGATACTGGTCTGGGTGGAGTCGCTTCTGAACATTGCCCTGTGCTTGTCTGA
 ATTGTACATGAAAAGGACTGGAGCAAGAAGGAAGTCCCAGGAATGGCAATGGAGTCACCTTAGAGCCA
 AGTGAAGCCAACATTAAACATGAGCGA

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR208965 protein sequence
 Red=Cloning site Green=Tags(s)


MGSTLGCHRSIPRPSDLSHNRKFSACNFSNILVNQERLNINTATEEELMTLPGVTRAVARSIVEYREY
 IGGFKKVEDLALVSGVGATKLEQVKFEICVSSKGNASQHSPLRRDLLAEQQPHHLTTTVPLTPRVNIN
 TATLAQLMSVRGLSEKMALSIDYRREHGPFSSVEDLVRMDGINAAFLDRIRHQVFAERSRPPSTHTNGG
 LTFTAKPHPSPTSLSLQSEDLDPGGPTQIISMRSVEAFGGMRDGRPVFRLATWNLQGCSEKANNPG
 VREVVCMTLLENSIKLLAVQELLDKEALEKFCTELNQPIPNIRKWKSGRCWRSIVAEPKSNQLQKGPC
 YSGFLWDTAANVELRDIPGRESSPSNGHAKAVGPSFPLTRFKVGSNDLTLVNLQLTALALPGAENSSKNH
 SDGHRLLNFALTLQETLKGEKDVVILGDFGQGPSNDYDILRREKFHHLVPAHTFTNISTRNPQGSKSVD
 NIWISLSLKKVFTGHWAVVREGLTNPWIPDNWSWGGVASEHCPVLAELYMEKDWSSKEVPRNGNGVTLEP
 SEANIKHER

TRTRPLE**QKLISEEDLAANDILDYKDDDDKV**

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_026189

ORF Size: 1707 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_026189.2 , NM_026189.3 , NP_080465.3
RefSeq Size:	2689 bp
RefSeq ORF:	1710 bp
Locus ID:	67484
UniProt ID:	Q3TGW2
Cytogenetics:	9 A4
MW:	63 kDa

Product images:



Circular map for MR208965