

Product datasheet for **MR208980L4V**

Haus3 (NM_146159) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Haus3
Synonyms:	D4S43h; D5H4S43; D5H4S43E
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_146159
ORF Size:	1422 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR208980).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_146159.1
RefSeq Size:	1980 bp
RefSeq ORF:	1713 bp
Locus ID:	231123
UniProt ID:	Q8QZX2
Cytogenetics:	5 17.84 cM



Gene Summary:

Contributes to mitotic spindle assembly, maintenance of centrosome integrity and completion of cytokinesis as part of the HAUS augmin-like complex.[UniProtKB/Swiss-Prot Function]