

## Product datasheet for **MR208995L3V**

### Zfp354a (BC087540) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Zfp354a (BC087540) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zfp354a
Synonyms:	AW488485; kid1; Tcf17; Znf354a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC087540
ORF Size:	2835 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208995).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC087540.1</a>
RefSeq Size:	4180 bp
RefSeq ORF:	1715 bp
Locus ID:	21408
Cytogenetics:	11 30.97 cM
Gene Summary:	It may play a role in renal development and may also be involved in the repair of the kidney after ischemia-reperfusion or folic acid administration.[UniProtKB/Swiss-Prot Function]



[View online »](#)