

Product datasheet for **MR208999L4V**

Ddx27 (BC011321) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ddx27 (BC011321) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ddx27
Synonyms:	C86129
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC011321
ORF Size:	1419 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR208999).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC011321 , AAH11321
RefSeq Size:	2483 bp
RefSeq ORF:	1718 bp
Locus ID:	228889
Cytogenetics:	2 H3
Gene Summary:	Probable ATP-dependent RNA helicase. Component of the nucleolar ribosomal RNA (rRNA) processing machinery that regulates 3' end formation of ribosomal 47S rRNA. [UniProtKB/Swiss-Prot Function]



[View online »](#)