

Product datasheet for **MR209013L3V**

Muc13 (NM_010739) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Muc13 (NM_010739) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Muc13
Synonyms:	14/A10; 114/A10; AI159736; Lrrp; Ly64; NJ-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010739
ORF Size:	1722 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209013).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010739.1 , NP_034869.1
RefSeq Size:	2969 bp
RefSeq ORF:	1722 bp
Locus ID:	17063
UniProt ID:	P19467
Cytogenetics:	16 B3
Gene Summary:	Epithelial and hemopoietic transmembrane mucin that may play a role in cell signaling. [UniProtKB/Swiss-Prot Function]



[View online »](#)