

Product datasheet for **MR209086L4V**

Flcn (NM_146018) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Flcn (NM_146018) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Flcn
Synonyms:	AU014660; B430214A04Rik; Bhd; FLCL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_146018
ORF Size:	1740 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209086).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_146018.1
RefSeq Size:	3460 bp
RefSeq ORF:	1740 bp
Locus ID:	216805
UniProt ID:	Q8QZS3
Cytogenetics:	11 B1.3
Gene Summary:	May be a tumor suppressor. May be involved in energy and/or nutrient sensing through the AMPK and mTOR signaling pathways. May regulate phosphorylation of RPS6KB1. [UniProtKB/Swiss-Prot Function]



[View online »](#)