

Product datasheet for **MR209101L4V**

Sppl2c (NM_001082535) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sppl2c (NM_001082535) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sppl2c
Synonyms:	4933407P14Rik; IMP5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001082535
ORF Size:	1743 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209101).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001082535.1 , NP_001076004.1
RefSeq Size:	2250 bp
RefSeq ORF:	1746 bp
Locus ID:	237958
UniProt ID:	A2A6C4
Cytogenetics:	11 E1
Gene Summary:	Intramembrane-cleaving aspartic protease (I-CLiP) that may be able to cleave type II membrane signal peptides in the hydrophobic plane of the membrane.[UniProtKB/Swiss-Prot Function]



[View online »](#)