

Product datasheet for **MR209113**

Ubqln1 (NM_026842) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ubqln1 (NM_026842) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ubqln1
Synonyms:	1110046H03Rik; 1810030E05Rik; AU019746; C77538; D13Erttd372e; Da41; Dsk2; Plic-1; Plic1; Xdrp1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR209113 representing NM_026842
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCGAGAGCGCAGAGAGCGGGCCACCGGGCGCGCAGGACAGTGGCGCCACGGCGCCCGCAG
 AGCCCAAGATCATGAAGTACACGGTGAAGACGCCTAAAGAGAAGGAGGAGTTGCCGTGCCGAGAACAG
 CTCGGTCCAGCAGTTCAAAGAGGAAATCTCAAACGTTTCAAATCTCATATTGACCAACTTGTGTTGATA
 TTTGCTGGAAAAATTTAAAAGATCAAGATACTCTGAGTCAGCATGGGATTCATGATGGACTTACGGTTC
 ACCTTGTGCATCAAACACAGAACAGGCCCAAGATAATTCAGCTCAGCAAACAAATGCCCTGGAAGCAC
 TGTGACCAGCTCACCAGCTCCCGATAGCAACCCACGCTGGTCTGCTGCTAACAGCTCCTTTGGTGT
 GGCGGACTTGGAGGACTGCAGGCTTAGCAGCTGGGTTTAAACACCACCAACTCTCTGAACTCCAGA
 GCCAAATGCAGCGGCAACTTTTGTCTAACCTGAAATGATGGTCCAGATCATGAAAAATCCCTTTGTCCA
 GAGCATGCTCTAAATCCTGACCTGATGCGGCAATTGATTATGGCCAATCCACAGATGCAACAGTTGATA
 CAGAGAAATCCAGAGATTAGTCATATGCTTAATAATCCAGATATAATGAGACAAACACTGGAAGTTGCGA
 GAAATCCAGCGATGATGCAGGAAATGATGAGAAACCAGGACCGAGCCTTGAGCAACCTAGAAAGTATTCC
 TGGGGGTATAATGCCTTACGGCGCATGTACCCGATATCCAGGAGCCTATGCTGAATGCTGCGCAAGAA
 CAGTTTGGTGGTAATCCATTTGCGTCTTAGTGAGCAGTTCGTCCTCAGCAGAAGGGACTCAGCCTTCTC
 GAACAGAAAAATAGGGATCCGCTCCCAACCCGTTGGGCTCCGCAGACTCCAGAGCTCACCAGCTCCGG
 TACCACTGGCAGCACCCTAACACCATGAGCACCTCTGGTGGCACTGCCACCAGTACCCCTGCTGGACAG
 AGCACCTCGGGCCCAAGCTTGGTGCCTGGCGCAGGAGCTAGTATGTTCAACACCCCGGAATGCAGAGCC
 TGCTCGCAGAGATAACTGAGAACCACAGCTCATGCAGAACATGCTGTCGGCCCATACATGAGGAGCAT
 GCTGCAGTCCCTGAGCCAGAACCCGACCTTGGTGCAGATGATGCTGAATAATCCCTATTTGCTGGA
 AATCCTCAGTTCAGAACAATGAGACACAGCTCCCAACTTCTCCAAACAATGCAGAACCCCGACA
 CGCTCTCAGCGATGTCGAACCTAGAGCCATGCAGGCGCTGCTGCAGATCCAGCAGGGCTTGCAGACGCT
 GGCCACAGAAGCTCCCGCCTCATCCAGGGTTTACTCCTGGCTTGGCAGCAGGCAATTCTGGAGGCTCT
 TCCGGAACCAATGCACCTAGCACTGCACCTAGTGAAGATACGAACCCCAAGGGGCACTGCTGAGCCAG
 GCCACCAGCAGTTTATCCAGCAGATGCTGCAGGCCCTTGGGGGGTGAATCCTCAGCTGCAGAGTCCAGA
 AGTCAGATTTAGCAACAACCTGGAACAACCTCAGTCCATGGGATCTTGAACCGTGAAGCAACCTGCAA
 GCTCTGATAGCCACAGGGGGCGACATCAATGCAGCAATCGAAAGTTGCTGGGCTCCAGCCGTCA

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR209113 representing NM_026842
 Red=Cloning site Green=Tags(s)

MAESAESGGPPGAQDSAADGGPAEPKIMKVTVKTPKEKEEFVAVPENSSVQQFKEEISKRFKSHIDQLVLI
 FAGKILKDQDTLSQHGIDHGLTVHLVIKQNRPDNSAQQTNPAGSTVTSSPAPDSNPTSGSAANSSFGV
 GGLGGLAGLSSLGLNTTNFSELQSQMQRQLLSNPMMVQIMENPFVQSMLSNPDLMRQLIMANPQMQLI
 QRNPEISHMLNPNPDIRQTLELARNPAMMQEMMRNQDRALSNLESIPGGYNALRRMYTDIQEPLNAAQE
 QFGGNPFASLVSSSSAEGTQPSRTENRDLPLPNWAPQTSQSSPASGTTGSTNTMSTSGGTATSTPAGQ
 STSGPSLVPGAGASMFNTPGMQSLLLQITENPQLMQNMLSAPYMRSMQLSLSQNPDLAAQMMLNPNPLFAG
 NPQLQEQRQLPTFLQQMQNPDTLSAMSNPRAMQALLQIQQLQLATEAPGLIPGFTPGLAAGNSGGG
 SGTNAPSTAPSEDTPNQQGTAEPGHQQFIQQMLQALAGVNPQLQSPEVRFQQQLQLSAGFLNREANLQ
 ALIATGGDINAAIERLLGSQPS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_026842

ORF Size: 1746 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026842.2](#)

RefSeq Size: 3693 bp

RefSeq ORF: 1749 bp

Locus ID: 56085


UniProt ID: [Q8R317](#)

Cytogenetics: 13 30.95 cM

MW: 62.4 kDa

Gene Summary: Plays an important role in the regulation of different protein degradation mechanisms and pathways including ubiquitin-proteasome system (UPS), autophagy and endoplasmic reticulum-associated protein degradation (ERAD) pathway. Mediates the proteasomal targeting of misfolded or accumulated proteins for degradation by binding (via UBA domain) to their polyubiquitin chains and by interacting (via ubiquitin-like domain) with the subunits of the proteasome. Plays a role in the ERAD pathway via its interaction with ER-localized proteins UBXLN4, VCP and HERPUD1 and may form a link between the polyubiquitinated ERAD substrates and the proteasome. Plays a role in unfolded protein response (UPR) by attenuating the induction of UPR-inducible genes, DDTI3/CHOP, HSPA5 and PDIA2 during ER stress. Involved in the regulation of macroautophagy and autophagosome formation; required for maturation of autophagy-related protein LC3 from the cytosolic form LC3-I to the membrane-bound form LC3-II and may assist in the maturation of autophagosomes to autolysosomes by mediating autophagosome-lysosome fusion. Negatively regulates the TICAM1/TRIF-dependent toll-like receptor signaling pathway by decreasing the abundance of TICAM1 via the autophagic pathway. Plays a key role in the regulation of the levels of PSEN1 by targeting its accumulation to aggresomes which may then be removed from cells by autophagocytosis. Promotes the ubiquitination and lysosomal degradation of ORA11, consequently downregulating the ORA11-mediated Ca²⁺ mobilization. Suppresses the maturation and proteasomal degradation of amyloid beta A4 protein (A4) by stimulating the lysine 63 (K63)-linked polyubiquitination. Delays the maturation of A4 by sequestering it in the Golgi apparatus and preventing its transport to the cell surface for subsequent processing (By similarity). Links CD47 to the cytoskeleton (PubMed:10549293).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR209113