

Product datasheet for **MR209117**

Ddi2 (BC031491) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ddi2 (BC031491) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ddi2
Synonyms:	1110056G13Rik; 1700027M01Rik; 9130022E05Rik; AI604911; AU040698
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR209117 representing BC031491
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCATCATTGCCGACTTCAGATGGGTTTGACCATCCAGCTCCTTCAGGGCAGAGTCTGAGGTTGGTA
 GCCCGACGAGTCTCGCTCGCTCTGTTTCTGCTTCGCTCTGCGCCATCAAGCCCGGTGACCCCAATAGCAT
 TGAATCTCTAGCTATGGAGGCTACGAAGCTTCAGCTGAATTCAGACAACTCTAAGAAAAACAGACCCCT
 CCTCTCTGCAGGTTCTTCTGACCTTGCTTCCTCAGCAGAGCAGAGTCTAGCCATGCCTTTCCATAAGT
 CATCAAAAGAAGCCGTTGTTGCAGGTAATCTGGAGAAATCTGTTGAGAAAGGAACCCAGGGCCTCAGAGT
 GTATCTCCACACAAGACAGGACGCTAGTTAACTCTCACAACACTGGGATGCGGAGCCACAGATATTT
 GCGGAGGAAAAGAGTTGGCATCCAGAAAATCAGACCCCAAGTCCCGTGAACGGCCTTCAGCAGCAGAG
 AAACAGGGAGTGTACAGCGAGAGGCTGGACAGCAGAGTGTCCACAGGACCAGGGCTGTCTTTGTACGC
 AGAAGACCTTGAGCTTCATGAAGAAGTTGTGAGTTTGAAGCTCTGAGGAAAGGTGAGCTACAAAGACAC
 GCTCATCTTCCCAGTGCAGAGAAGGGTCTTCCAGCTTCAGGACTCTGTAGCTGTCCATGCTCAGAAGCCC
 TGATGGAAGTAGATACAGCTGAACAGTCTCTGGTTGCTATGTGCAGCTCAACAGGCAGGCAGGATGCCGT
 CATCAAGAGCCCTAGTGTAGCACATCTCGCTTCAGATAATCCCACTATGGAAGTGAGACATTACAGTCT
 AACCCGTCCTGTGAGCCTGTGGAACATTCATCTTGACTCGGGAATTGCAGCTCCAGAAGATAATGTTG
 ACATGTCTACAATGGATAACAAAGATGACAATTCCTCTTCCCTTCTAAGTGGCCACGGTCAGCCCTCTGT
 GGAATCAGCAGAAGAATTTTGTTCATCTGTACAGTGGCCTTGAAAGAACTCCATGAGCTTTTGGTCATT
 AGCTGTAAGCCAGCTTCAGAAGAGTCACTGAGCATGTTACCTGTCAGTCAGAGATAGGAGCTGAGAGCC
 AACCAAGTGTTCAGACCTTTTCAAGGAAGAGGTCCAAAGTGTGATTTGACCCCTAGTGACCACTACTC
 ACAAGGCTCTGTGACAGGCCACCTCTGAATCAGGAAAGACAGAAATCGTAGGAAGTGGCCCTTGTGCT
 GCGGTAGAAGATGAAGCATCCACTAGCTTTGAAGGTCTGGGTGATGGCTGTACCTGACCGAGAAGATG
 TCCGCAGATCAACAGAGTCAGCCAGGAAGAGCTGTTCTGTGCGCCATAACCTCAGCTAACTGTCTGAGCA
 GTTGCCCTGCACCTCAGGGGTAGAAATAGCACCTGAACTTGACAGCAAGTGAAGGTGCCACAGTCAGCCT
 TCAGAGCATGTGCATAATCCAGGCCCAGACAGGCCGAGACCAGCAGTGTCTGCCCTGGAGCTGGGTTC
 CCCGTAGTGGATTGGACCAACCTCCACACAGTCTTGTCCACCCCTCCGTTCTTCCACCGTTCATCTT
 TCCTGCTGCAGATGTTGACAGGATCTTGGTGCCGGCTTCACTCTGCAGGAAGCGCTCGGGGCTCTGCAT
 CGAGTTGGTGGGAATGCAGACCTTGCCCTTCTGTTTTGTAGCAAAGAACATTGTAGTCCCAACA

ACGCGTACGCGGCGGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR209117 representing BC031491
 Red=Cloning site Green=Tags(s)


MSSLPTSDGFDHPAPSGQSPEVGSPTSLARSVSASVCAIKPGDPNSIESLAMEATKASAEFQTNSSKTD
 PPLQVLPDLASSAEQSLAMPFHKSSKEAVVAGNLEKSVEKGTQGLRVYLHTRQDASLTLTTGMREPQIF
 AEEKSWHPENQTPSPVNLQQHRETGSVQREAGQSVQDQGCCLDAEDLELHEEVVSLEALRKELQRH
 AHLPSAEKGLPASGLCSCPCSEALMEVDTAESLAVAMCSSTGRQDAVIKSPVAHLASDNPTMEVETLQS
 NPSCPEVHSILTRELQLPEDNVDMDNDKDDNSSLLSGHGQPSVESAEFCSSVTVALKELHELLVI
 SCKPASEESPEHVTCQSEIGAESQPSVSDLSGRRVQSVHLTPSDQYSQGSCHQATSESGKTEIVGTAPCA
 AVEDEASTSFEGLGDGLSPDREDVRRSTESARKSCSVAITSAKLSEQLPCTSGVEIAPELAASEGAHSQP
 SEHVHNPGPDRPETSSVCPGAGLPRSGLDQPPTQSLSTPSVLPPIFPAAADVDRILGAGFTLQEALGALH
 RVGGNADLALLVLLAKNIVVPT

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: BC031491

ORF Size: 1746 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC031491](#)


RefSeq Size: 2307 bp

RefSeq ORF: 1748 bp

Locus ID: 68817

Cytogenetics:	4 E1
MW:	84.6 kDa
Gene Summary:	Aspartic protease that mediates the cleavage of NFE2L1/NRF1 at 'Leu-104', thereby promoting release of NFE2L1/NRF1 from the endoplasmic reticulum membrane. Ubiquitination of NFE2L1/NRF1 is a prerequisite for cleavage, suggesting that DDI2 specifically recognizes and binds ubiquitinated NFE2L1/NRF1. Seems to act as a proteasomal shuttle which links the proteasome and replication fork proteins like RTF2. Required, with DDI1, for cellular survival following replication stress. Together or redundantly with DDI1, removes RTF2 from stalled forks to allow cell cycle progression after replication stress and maintains genome integrity. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR209117