

Product datasheet for **MR209118L4V**

Pxk (NM_145458) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pxk (NM_145458) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pxk
Synonyms:	C230080L11Rik; D14Ertd813e; MONaKA; Monaka; Slob
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145458
ORF Size:	1749 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209118).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145458.3
RefSeq Size:	2807 bp
RefSeq ORF:	1749 bp
Locus ID:	218699
UniProt ID:	Q8BX57
Cytogenetics:	14 4.7 cM
Gene Summary:	Binds to and modulates brain Na,K-ATPase subunits ATP1B1 and ATP1B3 and may thereby participate in the regulation of electrical excitability and synaptic transmission. May not display kinase activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)