

Product datasheet for **MR209122L3V**

Asb15 (NM_080847) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Asb15
Synonyms:	4930400E23Rik; Asb-15
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_080847
ORF Size:	1749 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR209122).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_080847.1 , NP_543123.2
RefSeq Size:	4268 bp
RefSeq ORF:	1752 bp
Locus ID:	78910
UniProt ID:	Q8VHS6
Cytogenetics:	6 A3.1



Gene Summary:

May be a substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.[UniProtKB/Swiss-Prot Function]