

## Product datasheet for **MR209175L3V**

### Cdk5rap1 (NM\_025876) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdk5rap1 (NM_025876) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cdk5rap1
Synonyms:	2310066P17Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025876
ORF Size:	1767 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209175).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_025876.2</a> , <a href="#">NP_080152.1</a>
RefSeq Size:	1970 bp
RefSeq ORF:	1767 bp
Locus ID:	66971
UniProt ID:	<a href="#">Q8BTW8</a>
Cytogenetics:	2 H1
Gene Summary:	Probable regulator of CDK5 activity. May inhibit CDK5 function via its interaction with CDK5R1 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)