

Product datasheet for **MR209188L3V**

Fam161b (NM_172581) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Fam161b
Synonyms:	9330169D17; 9830169C18Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172581
ORF Size:	1767 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR209188).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_172581.2
RefSeq Size:	3357 bp
RefSeq ORF:	1770 bp
Locus ID:	217705
UniProt ID:	Q8CB59
Cytogenetics:	12 D1

