

Product datasheet for **MR209194L3V**

Rangap1 (NM_011241) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rangap1 (NM_011241) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rangap1
Synonyms:	C79654; Fug1; mKIAA1835
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011241
ORF Size:	1770 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209194).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011241.3
RefSeq Size:	3002 bp
RefSeq ORF:	1770 bp
Locus ID:	19387
UniProt ID:	P46061
Cytogenetics:	15 38.21 cM



[View online »](#)

Gene Summary:

GTPase activator for RAN. Converts cytoplasmic GTP-bound RAN to GDP-bound RAN, which is essential for RAN-mediated nuclear import and export (PubMed:18305100). Mediates dissociation of cargo from nuclear export complexes containing XPO1, RAN and RANBP2 after nuclear export (By similarity). Required for postimplantation embryonic development (PubMed:8314081).[UniProtKB/Swiss-Prot Function]