

Product datasheet for MR209200L4

Ifnar1 (NM_010508) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ifnar1 (NM_010508) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Ifnar1
Synonyms:	CD118; Ifar; Ifnar; Ifrc; Infar
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209200).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_010508
ORF Size:	1770 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_010508.1 , NP_034638.1
RefSeq Size:	3894 bp
RefSeq ORF:	1773 bp
Locus ID:	15975
UniProt ID:	P33896
Cytogenetics:	16 52.98 cM
Gene Summary:	Component of the receptor for type I interferons, including interferons alpha, IFNB1 and IFNW1 (PubMed:1533935, PubMed:14532120, PubMed:23872679). Functions in general as heterodimer with IFNAR2 (By similarity). Type I interferon binding activates the JAK-STAT signaling cascade, and triggers tyrosine phosphorylation of a number of proteins including JAKs, TYK2, STAT proteins and the IFNR alpha- and beta-subunits themselves (PubMed:14532120). Can form an active IFNB1 receptor by itself and activate a signaling cascade that does not involve activation of the JAK-STAT pathway (PubMed:23872679). Contributes to modulate the innate immune response to bacterial lipopolysaccharide (PubMed:23872679).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR209200L4