

Product datasheet for **MR209215L3V**

Lrrtm4 (BC037216) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lrrtm4 (BC037216) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lrrtm4
Synonyms:	7530419J18Rik; A230052N11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC037216
ORF Size:	1773 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209215).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC037216.1
RefSeq Size:	2976 bp
RefSeq ORF:	1775 bp
Locus ID:	243499
Cytogenetics:	6 C3
Gene Summary:	May play a role in the development and maintenance of the vertebrate nervous system. Exhibits strong synaptogenic activity, restricted to excitatory presynaptic differentiation (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)