

Product datasheet for **MR209244L4V**

HnrpII (BC012849) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HnrpII (BC012849) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	HnrpII
Synonyms:	PHGPx
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC012849
ORF Size:	1782 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209244).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC012849 , AAH12849
RefSeq Size:	3017 bp
RefSeq ORF:	1784 bp
Locus ID:	72692
Cytogenetics:	17 E3
Gene Summary:	RNA-binding protein that functions as regulator of alternative splicing for multiple target mRNAs, including PTPRC/CD45 and STAT5A. Required for alternative splicing of PTPRC. [UniProtKB/Swiss-Prot Function]


[View online »](#)