

Product datasheet for **MR209250L4V**

Setd3 (BC019973) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Setd3 (BC019973) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Setd3
Synonyms:	2610102I01Rik; 2610305M23Rik; 6030490B17Rik; D12ErtD771e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC019973
ORF Size:	1782 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209250).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC019973.1
RefSeq Size:	2477 bp
RefSeq ORF:	1784 bp
Locus ID:	52690
Cytogenetics:	12 59.23 cM


[View online »](#)

Gene Summary:

Protein-histidine N-methyltransferase that specifically mediates methylation of actin at 'His-73' (PubMed:30626964). Histidine methylation of actin is required for smooth muscle contraction of the laboring uterus during delivery (PubMed:30626964). Does not have protein-lysine N-methyltransferase activity and probably only catalyzes histidine methylation of actin (PubMed:30626964).[UniProtKB/Swiss-Prot Function]