

Product datasheet for **MR209261L4V**

6330514A18Rik (BC099679) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	6330514A18Rik (BC099679) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	6330514A18Rik
Synonyms:	6330514A18Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC099679
ORF Size:	1785 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209261).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC099679 , AAH99679
RefSeq Size:	2089 bp
RefSeq ORF:	1787 bp
Locus ID:	216166
Cytogenetics:	10 C1
Gene Summary:	Inactive serine/threonine-protein kinase that plays a role in cell cycle progression and neuronal differentiation.[UniProtKB/Swiss-Prot Function]



[View online »](#)