

## Product datasheet for **MR209264L3V**

### **D15Erttd621e (BC033609) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	D15Erttd621e (BC033609) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	D15Erttd621e
Synonyms:	mKIAA0493, MGC40987
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC033609
ORF Size:	1785 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209264).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC033609</a> , <a href="#">AAH33609</a>
RefSeq Size:	4689 bp
RefSeq ORF:	1787 bp
Locus ID:	210998
Cytogenetics:	15 24.76 cM
Gene Summary:	As component of the WDR11 complex acts together with TBC1D23 to facilitate the golgin-mediated capture of vesicles generated using AP-1.[UniProtKB/Swiss-Prot Function]


[View online »](#)