

## Product datasheet for **MR209289L4V**

### Zfp143 (BC037658) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Zfp143 (BC037658) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zfp143
Synonyms:	AA959806; AU015869; D7Erttd805e; KRAB14; pHZ-1; SBF; Staf; Zfp79; Zfp80-rs1; Znf143
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC037658
ORF Size:	1794 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209289).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC037658.1</a>
RefSeq Size:	2980 bp
RefSeq ORF:	1796 bp
Locus ID:	20841
Cytogenetics:	7 57.7 cM
Gene Summary:	Transcriptional activator. Activates the gene for selenocysteine tRNA (tRNA <sup>sec</sup> ). Binds to the SPH motif of small nuclear RNA (snRNA) gene promoters. Participates in efficient U6 RNA polymerase III transcription via its interaction with CHD8 (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)