

Product datasheet for **MR209300L3V**

Cel (NM_009885) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cel (NM_009885) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cel
Synonyms:	1810036E18Rik; BAL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009885
ORF Size:	1800 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209300).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009885.1 , NP_034015.1
RefSeq Size:	2046 bp
RefSeq ORF:	1800 bp
Locus ID:	12613
UniProt ID:	Q64285
Cytogenetics:	2 19.38 cM
Gene Summary:	Catalyzes fat and vitamin absorption. Acts in concert with pancreatic lipase and colipase for the complete digestion of dietary triglycerides (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)