

## Product datasheet for **MR209305L3V**

### Slc13a3 (NM\_054055) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc13a3 (NM_054055) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc13a3
Synonyms:	NaDC-3; NaDC3; SDCT2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_054055
ORF Size:	1803 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209305).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_054055.2</a> , <a href="#">NP_473396.1</a>
RefSeq Size:	3224 bp
RefSeq ORF:	1803 bp
Locus ID:	114644
UniProt ID:	<a href="#">Q91Y63</a>
Cytogenetics:	2 H3
Gene Summary:	High-affinity sodium-dicarboxylate cotransporter that accepts a range of substrates with 4-5 carbon atoms. The stoichiometry is probably 3 Na(+) for 1 divalent succinate. [UniProtKB/Swiss-Prot Function]



[View online »](#)