

Product datasheet for **MR209308L3V**

Pola2 (NM_008893) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pola2 (NM_008893) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pola2
Synonyms:	AI573378
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008893
ORF Size:	1803 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209308).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008893.2 , NP_032919.2
RefSeq Size:	2424 bp
RefSeq ORF:	1803 bp
Locus ID:	18969
UniProt ID:	P33611
Cytogenetics:	19 A
Gene Summary:	May play an essential role at the early stage of chromosomal DNA replication by coupling the polymerase alpha/primase complex to the cellular replication machinery.[UniProtKB/Swiss-Prot Function]



[View online »](#)