

Product datasheet for **MR209322L4V**

Ptgs1 (NM_008969) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ptgs1 (NM_008969) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ptgs1
Synonyms:	COX; Cox-; Cox-1; Cox-3; COX1; Pgh; PGHS-1; Pghs1; PHS 1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008969
ORF Size:	1806 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209322).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008969.3
RefSeq Size:	2864 bp
RefSeq ORF:	1809 bp
Locus ID:	19224
UniProt ID:	P22437
Cytogenetics:	2 24.19 cM



[View online »](#)

Gene Summary:

This is one of two genes encoding similar enzymes that catalyze the conversion of arachinodate to prostaglandin. The encoded protein regulates angiogenesis in endothelial cells, and is inhibited by nonsteroidal anti-inflammatory drugs such as aspirin. Based on its ability to function as both a cyclooxygenase and as a peroxidase, the encoded protein has been identified as a moonlighting protein. [provided by RefSeq, Jan 2014]