

## Product datasheet for **MR209330L4V**

### Galnt10 (NM\_134189) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Galnt10 (NM_134189) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Galnt10
Synonyms:	AU018154; C330012K04Rik; GalNAc-T9; GalNAc-T10; Galnt9; ppGaNTase
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_134189
ORF Size:	1812 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209330).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_134189.2</a> , <a href="#">NP_598950.2</a>
RefSeq Size:	4725 bp
RefSeq ORF:	1812 bp
Locus ID:	171212
UniProt ID:	<a href="#">Q6P9S7</a>
Cytogenetics:	11 B1.3



[View online »](#)

**Gene Summary:**

Catalyzes the initial reaction in O-linked oligosaccharide biosynthesis, the transfer of an N-acetyl-D-galactosamine residue to a serine or threonine residue on the protein receptor. Has activity toward Muc5Ac and EA2 peptide substrates (By similarity).[UniProtKB/Swiss-Prot Function]