

Product datasheet for **MR209334L3V**

Gga2 (NM_028758) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gga2 (NM_028758) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gga2
Synonyms:	1200007E24Rik; mKIAA1080
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_028758
ORF Size:	1812 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209334).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_028758.2 , NP_083034.1
RefSeq Size:	4891 bp
RefSeq ORF:	1812 bp
Locus ID:	74105
UniProt ID:	Q6P5E6
Cytogenetics:	7 F2



[View online »](#)

Gene Summary:

Plays a role in protein sorting and trafficking between the trans-Golgi network (TGN) and endosomes. Mediates the ARF-dependent recruitment of clathrin to the TGN and binds ubiquitinated proteins and membrane cargo molecules with a cytosolic acidic cluster-dileucine (DXLL) motif. Mediates export of the GPCR receptor ADRA2B to the cell surface. Regulates retrograde transport of phosphorylated form of BACE1 from endosomes to the trans-Golgi network.[UniProtKB/Swiss-Prot Function]