

## Product datasheet for **MR209347L4V**

### **Ncam1 (BC011310) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ncam1 (BC011310) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ncam1
Synonyms:	Ncam, CD56, NCAM-1, NCAM-120, NCAM-140, NCAM-180, E-NCAM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC011310
ORF Size:	1815 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209347).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC011310.1</a>
RefSeq Size:	2544 bp
RefSeq ORF:	1817 bp
Locus ID:	17967
Cytogenetics:	9 26.83 cM
Gene Summary:	This protein is a cell adhesion molecule involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc.[UniProtKB/Swiss-Prot Function]



[View online »](#)