

Product datasheet for **MR209399L4V**

Tomm70a (NM_138599) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tomm70a (NM_138599) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tomm70a
Synonyms:	2610044B22Rik; D16lum22; D16lum22e; D16Wsu109e; mKIAA0719; Tom70; Tomm70
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138599
ORF Size:	1833 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209399).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138599.2 , NP_613065.2
RefSeq Size:	3778 bp
RefSeq ORF:	1836 bp
Locus ID:	28185
UniProt ID:	Q9CZW5
Cytogenetics:	16 34.22 cM
Gene Summary:	Receptor that accelerates the import of all mitochondrial precursor proteins. [UniProtKB/Swiss-Prot Function]



[View online »](#)