

## Product datasheet for **MR209404L4V**

### Igsf8 (NM\_080419) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Igsf8 (NM_080419) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Igsf8
Synonyms:	ESTM34; EWI-2; KCT-4; PGRL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080419
ORF Size:	1836 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209404).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_080419.1</a>
RefSeq Size:	2256 bp
RefSeq ORF:	1836 bp
Locus ID:	140559
UniProt ID:	<a href="#">Q8R366</a>
Cytogenetics:	1 79.59 cM



[View online »](#)

**Gene Summary:**

May play a key role in diverse functions ascribed to CD81 and CD9 such as oocytes fertilization or hepatitis C virus function. May regulate proliferation and differentiation of keratinocytes. May be a negative regulator of cell motility: suppresses T-cell mobility coordinately with CD81, associates with CD82 to suppress prostate cancer cell migration, regulates epidermoid cell reaggregation and motility on laminin-5 with CD9 and CD81 as key linkers. May also play a role on integrin-dependent morphology and motility functions. May participate in the regulation of neurite outgrowth and maintenance of the neural network in the adult brain.[UniProtKB/Swiss-Prot Function]