

Product datasheet for **MR209405**

LOC626578 (BC114209) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LOC626578 (BC114209) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LOC626578
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR209405 representing BC114209
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGACCCAACCACAAATGGCTCCCATTTGCCTTGTGAAAACCACAATGAACAGCTGTCAGTGAACCAGG
 AAGCCATAGAGATTCTGGACAAGATTTCTCAGCCAGTGGTAGTCGTGGCTATTGTTGGATGGTCCCGTAC
 TGGGAAGTCTATTTGATGAAGTGTCTGGCGGGACAGAATCACGGCTTCCTCTGGGCTCCACTGTGCAG
 TCTCAAACCAAGGGCATCTGGATGTGGTGCATGCCACACCCACCAAGCCAGAGCACACCTGGTCTCC
 TGGACACCGAGGGCTGGGGATGTGAAAAGGGTGATCCTAAGAAGCACTTGTGGATCTTGGCCCTCAG
 TGTGCTTCTGAGCAGCCTTCATCTACAACAGCATGATCACCATCAACCACAGGCCCTGGAGCAGCTG
 CATTATGTCACAGAACTCACAGAGCTGATCAGAGCAAAGTCTTCCCAAATCCTGCTGGAATAAAGAATT
 CCACAGAGTTTGTGAGTTTCTTCCAGACTTTGTCTGGATTGTTGGGATTTTCATGCTGGAGCTGAAGTT
 AAATGGGAAGACATACAAGTGTGACTACCTGGAGAATGCCTGAAGCTGATCCAGGTGACAAACCC
 AGAATGCAAGCATCCAATTCATGCAGGGAATGCATCAGACTTTTCTTTCTAACCAGGAAGTGTGTTGTCT
 TTGACCGGCCAACGCATGACAAAGAACTTTTACAAAACTTGATTCTATCACAGAAGACCAACTGGATCC
 TAAATTCAGGAAGTAACAAAGGCTTTTGTCTTACATCTTCACTTATGCCAAGATCAAGACCTAAAA
 GAGGGAATTAAGGTCACTGGGAATAAAGTAGGGATTCTGGTGACAACCTATGTGGATGCCATCAACAGTG
 GAGCAGTGCCCTTGTCTGGATGATGCTGTGACAACCTCTGGCCAGCGTGAGAACTCAGTAGCTGTGCAGAA
 AGCAGCCGACCACTATAGTGAGCAGATGGCCAGCGACTGAGGCTTCTACAGAAACGCTCCAGGAGCTG
 CTGGATGTGCATGCAGCCTGCGAGAAGGAAGCCATGGCTGTCTTCATGGAGTATTCCTTCAAGGATGAAA
 ATCAGCAATTCCTGAAGAAGCTGGTGAAATTAATAGGAGAGAACAAGAGCTTTTCTGTGCAAGAATGA
 AGAGGCATCAAAATAAATACTGTCAAGAAGAACTGGATCGACTTTTCAAAGGATTTTATGGAAAAATTTTCA
 ACATTTTTTGTCTTGTGGACACAAGCTTTACATGGACAAGAGGGAGAAGATTGAGCATGACTACTGGC
 AGGTTCCAGGAAGGGGTGAAGGCAAGTGAAGTCTTCCAGAGCTTTCTGCAGTCACAGGCCTTCATCGA
 GAGTTCCATCTTGCAGGCAGATACAGCGCTCACTGCTGGGAAGAAGGCCATTGCAGAGGAGCGTGCCAG
 AAGGTGGCGGCAGAGAAGGAGCAGGATCTGCTAAGACAGAAGCAGAAGGAGCAGCAGGAGTATATGGAGG
 CTCAAGAGAAAAGTCACAAGGAAAACCTAGAGCAACTGAGGAGAAAGCTGATGCAGGAGAGAGAGCAGCT
 CATCAAAGACCATAACATGATGTGAAGAAGCAGCTGAAGGATCAAAAGGCTTTCTTGAAGGATTT
 AAGAAGAAAGCTGAGGAAATGAACAAAGAGATACAGCAACTGAGAGATGTCATCAAGGATAAGAAAAGAA
 ACACTGATCGAATTAAGGATGCTCTCTTAATGGATTTTCTACAGTTCTTTTTCATTACCTGTCCGTTA
 TCTAAAGCATTTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR209405 representing BC114209
 Red=Cloning site Green=Tags(s)

MTQPQMAPICLVENHNEQLSVNQEAIEILDKISQPVVVVAIVGWSRTGKSYLMNCLAGQNHFPLGSTVQ
 SQTKGIMWCMHPHTKPEHTLVLLDTEGLGDVEKGDPKNDLWIFALSVLLSSTFIYNSMITINHQAELQ
 HYVTELTELIRAKSSPNPAGIKNSTEFVSFFPDFVWIVRDFMLELKLNGEDITSDDYLENALKLIPGDKP
 RMQASNSCRECIRLFFPNRKCFFVDRPTHDKELLQKLDSITEDQLDPKFQEVTKAFVSYIFTYAKIKTLK
 EGIKVTGNKLGILVTTYVDAINS GAVPCLDDAVTTLAQRENSVAVQKAADHYSEQMAQRLRLPTETLQEL
 LDVHAACEKEAMAVFMEYSFKDENQQFLKKLVELIGENKELFLSKNEEASNKYCQEELDRLSKDFMENIS
 TFFVPCGHKL YMDKREKIEHDYWQVPRKGVKASEVFQSFLQSQAFIESSILQADTALTAGKKAIAEERAQ
 KVAAEKEQDLLRQKQKEQQEYMEAQEKSHKENLEQLRRKLMQEREQLIKDHNMLKKQLKDQKAFLEEGF
 KKKAEE MNKEIQQLRDVIKDKKRNTDRIKDALLNGFSTVL FHYLVRYL KHL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mm9073_b11.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: BC114209

ORF Size: 1833 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: BC114209.1

RefSeq Size: 2631 bp

RefSeq ORF: 1835 bp

Locus ID: 626578
 Cytogenetics: 5 E5
 MW: 70.4 kDa

Product images:



Circular map for MR209405