

## Product datasheet for **MR209458L4V**

### Adck2 (NM\_178873) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Adck2 (NM_178873) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Adck2
Synonyms:	AA408112; BB006057
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178873
ORF Size:	1854 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209458).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_178873.2</a> , <a href="#">NP_849204.1</a>
RefSeq Size:	3539 bp
RefSeq ORF:	1854 bp
Locus ID:	57869
UniProt ID:	<a href="#">Q6NSR3</a>
Cytogenetics:	6 B1
Gene Summary:	The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).[UniProtKB/Swiss-Prot Function]



[View online »](#)