

Product datasheet for **MR209492L4V**

Gbp3 (NM_018734) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gbp3 (NM_018734) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gbp3
Synonyms:	AW228655; GBP-3; GBP-4; Gbp4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018734
ORF Size:	1863 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209492).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018734.3
RefSeq Size:	2639 bp
RefSeq ORF:	1863 bp
Locus ID:	55932
UniProt ID:	Q61107
Cytogenetics:	3 H1
Gene Summary:	Binds GTP, GDP and GMP. Hydrolyzes GTP very efficiently; GDP rather than GMP is the major reaction product. Plays a role in erythroid differentiation.[UniProtKB/Swiss-Prot Function]



[View online »](#)