

Product datasheet for **MR209512L4V**

Ddx41 (NM_134059) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ddx41 (NM_134059) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ddx41
Synonyms:	2900024F02Rik; AA958953; ABS; AI324246
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_134059
ORF Size:	1869 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209512).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_134059.1 , NP_598820.1
RefSeq Size:	2213 bp
RefSeq ORF:	1869 bp
Locus ID:	72935
UniProt ID:	Q91VN6
Cytogenetics:	13 B1
Gene Summary:	Probable ATP-dependent RNA helicase. Is required during post-transcriptional gene expression. May be involved in pre-mRNA splicing.[UniProtKB/Swiss-Prot Function]



[View online »](#)