

Product datasheet for **MR209521L3V**

Xpnpep1 (NM_133216) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Xpnpep1 (NM_133216) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Xpnpep1
Synonyms:	D230045I08Rik; sAMP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133216
ORF Size:	1998 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR209521).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133216.3 , NP_573479.3
RefSeq Size:	2521 bp
RefSeq ORF:	2001 bp
Locus ID:	170750
UniProt ID:	Q6P1B1
Cytogenetics:	19 D2
Gene Summary:	Contributes to the degradation of bradykinin. Catalyzes the removal of a penultimate prolyl residue from the N-termini of peptides, such as Arg-Pro-Pro (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)